



# BLANKET WARMING CABINETS

FOR SAFE WARMING OF COTTON BLANKETS

# UK MANUFACTURED CABINETS

QED Scientific are proud to offer the BW-ECO blanket warming cabinet range,

Our blanket warming models are perfect for warming cotton blankets in operating theatre recovery, A&E or a maternity setting.

They offer instant warmth to patients in theatre recovery or who may be hypothermic in A&E.

With PID digital control, three levels of temperature safety and automatic overheat protection they offer the safest answer for warming cotton blankets and drapes.

Our Blanket Warming Cabinets can be combined with our ECO Fluid Warming cabinets to create a convenient space-saving combined unit.

## FEATURES

- 11 standard sizes of fluid warmer available, all having digital PID controllers
- Temperatures up to 80°C
- Automatic setting of the overheat safety cut-out.
- Secondary independent overheat safety fail-safe protection
- Capacities range from 4 up to about 50 cotton blankets in baskets
- Medical grade stainless steel interior with powder painted exterior
- Quiet internal fan ensures optimum temperature control with fast recovery after door opening
- Accessories include door locks, audible alarms, and mobile stands
- Custom-built to your specific size if required
- Combination cabinets for storing fluids and blankets are available
- CE Certified
- Manufactured in the UK
- 5 year UK warranty

Contact us to discuss your requirements and specification or to request further information.

### Office

[info@qedscientific.com](mailto:info@qedscientific.com)

01663 735494

# OUR RANGE

In addition to blanket warming, we also offer a wide range of standard sized cabinets for warming either fluids or contrast media as well as our bespoke sizes for unique projects.

## FLUID WARMERS

Designed to hold bottles from 1 litre up to 3 litres, or bags from half a litre to 3 litres.

Capacities from 9 to 240 x 1 litre bottles on heavy duty wire shelves; 2 to 60 x 3 litre bags in wire baskets. These have a temperature range from above ambient to 60°C.

Accessories include stands for for under bench models, lockable wheels for larger cabinets and door locks to prevent unauthorised access.

## COMBINATION WARMERS

For smaller operating theatres or A&E units where warm fluids and blankets are required but floor space is limited.

These are based on a combination of a fluid warmer stacked on top of a blanket warmer and fixed together.

Lockable wheels are fitted to all models.

Common accessories include digital door locks.

## CONTRAST MEDIA WARMERS

For warming bottles of contrast media in CT radiography rooms.

Models are available to warm from 10 to 120 x 50ml bottles stored on perforated shelves.

Bottles up to 200ml in size can be warmed and stored.

Door locks are often requested to be fitted to these cabinets.

Accessories also include wall fixing brackets and wall shelves.



