



FLUID WARMING CABINETS

SAFE, SIMPLE, EFFECTIVE

UK MANUFACTURED CABINETS

QED Scientific are proud to offer the FW-ECO fluid warming cabinets, perfect for warming bottles and bags of saline and other fluids in an operating theatre or maternity setting.

They are also known as solution warmers, and heated lotion cabinets.

With PID digital control, three levels of temperature safety and automatic overheat protection they offer the safest answer for warming fluids.

The chambers are medical grade polished stainless steel for easy cleaning to aid infection control procedures.

Our Blanket Warming Cabinets can be combined with our ECO Fluid Warming cabinets to create a convenient space-saving combined unit.

FEATURES

- 15 standard sizes of fluid warmer available, bespoke sizes for individual requirements
- Digital microprocessor control
- Temperatures up to 60°C, factory set to 37°C
- Automatic setting of the overheat safety cut-out. Secondary independent overheat safety fail-safe protection
- Capacities range from 6 half litre bottles up to 200 1-litre bottles or equivalent
- Medical grade stainless steel interior with powder painted exterior
- Quiet internal fan ensures optimum temperature control with fast recovery after door opening
- Accessories include baskets for bags of fluid, door locks, alarms, and stands
- Smaller models can be wall mounted
- Custom-built to your specific requirements if necessary
- Combined fluids and blanket cabinets are also available
- CE Certified
- Manufactured in the UK.
- 5 year UK warranty

Contact us to discuss your requirements and specification or to request further information.

Office

info@qedscientific.com

01663 735494

OUR RANGE

In addition to fluid warming, we also offer a wide range of standard sized cabinets for warming either blankets or contrast media as well as our bespoke sizes for unique projects.

BLANKET WARMERS

For warming cotton blankets in theatre recovery, A&E resus rooms and critical care areas. Cabinets are supplied with baskets for storing from 4 to 50 standard cotton blankets. The temperature for warming cotton blankets is limited to 80°C. Accessories include stands or wheels for smaller models, door locks to prevent unauthorised access

COMBINATION WARMERS

For smaller operating theatres or A&E units where warm fluids and blankets are required but floor space is limited.

These are based on a combination of a fluid warmer stacked on top of a blanket warmer And fixed together.

Lockable wheels are fitted to all models.

Common accessories include digital door locks.

CONTRAST MEDIA WARMERS

For warming bottles of contrast media in CT radiography rooms.

Models are available to warm from 10 to 120 x 50ml bottles stored on perforated shelves. Bottles up to 200ml in size can be warmed and stored.

Door locks are often requested to be fitted to these cabinets.

Accessories also include wall fixing brackets and wall shelves.



