



CONTRAST MEDIA WARMING CABINETS

FOR SAFE WARMING OF BOTTLES

UK MANUFACTURED CABINETS

QED Scientific are proud to offer the MW-ECO contrast media warming cabinet range,

Our MW-ECO models are designed specifically for warming bottles of contrast media in CT scanner suites.

They have a PID digital controller, three levels of temperature safety and automatic overheat protection.

They are factory set to control constantly at 37°C

There are a range of accessories/ options available depending on your needs..

As well as our standard sizes, we can manufacture a bespoke size if required.

FEATURES

- 7 standard sizes of media warmer are available, all with digital microprocessor control
- Temperatures up to 40°C, factory set to 37°C
- Automatic setting of the overheat safety cut-out. Secondary independent overheat safety fail-safe protection
- Capacities range from 10 to 140 bottles of media
- Medical grade stainless steel interior with powder painted exterior
- Quiet internal fan ensures optimum temperature control with fast recovery after door opening
- Accessories include an audible alarm and digital door lock
- Smaller models can be wall mounted
- Custom-built to your specific requirements if necessary
- CE Certified
- Manufactured in the UK.
- 5 year UK warranty

Contact us to discuss your requirements and specification or to request further information.

Office

info@qedscientific.com

01663 735494

OUR RANGE

In addition to contrast media, we also offer a wide range of standard sized cabinets for warming either fluids or cotton blankets as well as our bespoke sizes for unique projects.

FLUID WARMERS

Designed to hold bottles from 1 litre up to 3 litres, or bags from half a litre to 3 litres. Capacities from 9 to 240 x 1 litre bottles on heavy duty wire shelves; 2 to 60 x 3 litre bags in wire baskets. These have a temperature range from above ambient to 60°C. Accessories include stands for for under bench models, lockable wheels for larger cabinets and door locks to prevent unauthorised access.

BLANKET WARMERS

For warming cotton blankets in theatre recovery, A&E resus rooms and critical care areas. Cabinets are supplied with baskets for storing from 4 to 50 standard cotton blankets. The temperature for warming cotton blankets is limited to 80°C. Accessories include stands or wheels for smaller models, door locks to prevent unauthorised access

COMBINATION WARMERS

For smaller operating theatres or A&E units where warm fluids and blankets are required but floor space is limited. These are based on a combination of a fluid warmer stacked on top of a blanket warmer and fixed together. And fixed together. Lockable wheels are fitted to all models. Common accessories include digital door locks.



