



BESPOKE BUILT OVENS

In addition to our standard sizes of ovens, we can easily accommodate customer requests for one off designs.

Unique sizes up to 6000 litres capacity, position of controllers, higher temperature operation, side to side airflow and access with a forktruck for palletted drums of fluid are just some of the requests we have facilitated in the past.

Contact us to discuss your requirement.



GENERAL PURPOSE OVENS

Our standard range of general purpose ovens are available in 15 standard sizes from 10 litre capacity up to 355 litres and feature as standard an aluminium clad steel chamber, dial thermostats, independent safety cut out and a 3 year warranty.

Optional digital PID controllers and our attention to airflow design ensure that temperature stability meets Class 1 ($\pm 3^{\circ}\text{C}$) standards on most sizes with larger capacities meeting Class 2 ($\pm 6^{\circ}\text{C}$). another popular option is fan assistance which ensures a uniform temperature throughout the whole chamber.

Other options to consider are a higher operating temperature upgrade to 500°C , a full stainless steel chamber, programmable multi-step controllers and inner glass doors and repositioning of the controls to suit the position of the oven.

MODEL	CAPACITY (L)	NO OF SHELVES	EXTERNAL (H x W x D)	INTERNAL (H x W x D)
GP-OV-10	10	1	415 x 375 x 350	150 x 275 x 243
GP-OV-20	20	2	540 x 375 x 382	275 x 275 x 275
GP-OV-30	30	2	570 x 415 x 422	305 x 315 x 315
GP-OV-40	40	2	670 x 415 x 422	404 x 315 x 315
GP-OV-50	50	2	615 x 590 x 422	325 x 490 x 315
GP-OV-75	75	2	625 x 590 x 567	333 x 490 x 460
GP-OV-100	100	3	735 x 590 x 567	445 x 490 x 460
GP-OV-120	120	3	735 x 590 x 659	445 x 490 x 552
GP-OV-150	150	3	845 x 590 x 659	555 x 490 x 552
GP-OV-180	180	4	957 x 590 x 659	667 x 490 x 552
GP-OV-210	210	4	957 x 590 x 749	667 x 490 x 642
GP-OV-240	240	4	1178 x 590 x 659	888 x 490 x 552
GP-OV-280	280	4	941 x 750 x 769	651 x 650 x 662
GP-OV-350	350	4	1104 x 750 x 769	814 x 650 x 662
GP-OV-355	355	4	1620 x 590 x 659	1345 x 490 x 552



Contact us to discuss your requirements and specification or to request further information.

Office

info@qedscientific.com

01663 735494

